

EXHIBIT A

C-18 1.5%	S/F 30%	APC 10 5/24 1.2	PDGF 24 19
--------------	------------	--------------------------	------------------

Hep Seph / C-18 column of 1

see det 1st - 1st fraction test (was 1.5 ml)

3 Tubes of 500ul each dry vial

1 Tube resuspended in 40ul 5 mM HCl + 40ul sample buffer

ran 20ul on well #2

HUVE S/F media dialyzed against 1 N HAc then

0.1 N HAc for total of 24 hrs.

500ul Tubes dry vial. resuspended 1 Tube

in 20ul 5 mM HCl plus 20ul sample buffer

used 20ul in well #4

HUVE S/F media Hep Seph / APC: 10 2 PDGF


column of 40 ml → 4 ml in 2nd fraction

Took 500ul of 2nd fraction + dry vial.

resuspended 1 Tube in 40ul 5 mM HCl +

40ul sample buffer. Ran 20ul on gel well #6

PDGF is 24 ng of creatine kinase reagent.



PROGRAM # 7
 REGION A: LL-UL= 0- 19 LCR= 0 BKG= 0 % 2 SIGMA= .2
 REGION B: LL-UL= 2- 19 LCR= 0 BKG= 0 % 2 SIGMA= .2
 TIME= 1.00 QIP= SIS SCR= B/A K= 1.000

#	S#	TIME	CPMA/K	%DEV	CPMB/K	%DEV	QIP	FLAGS	SCR	MIN
1	1.00	30651	1.14	29698	1.16	16.0	32	.969	1	
2	1.00	13448	1.72	12950	1.76	15.9	32	.963	3	
3	1.00	15974	1.58	15408	1.61	15.8	32	.965	4	
4	1.00	11639	1.85	11190	1.89	15.9	32	.961	5	
5	1.00	3431	3.41	3249	3.51	14.8	32	.947	6	
6	1.00	10127	1.99	9795	2.02	16.1	27	.967	8	
7	1.00	60448	.81	58777	.82	16.4	32	.972	9	
8	1.00	39864	1.00	38734	1.02	16.4	27	.972	10	
9	1.00	45614	.94	44348	.95	16.8	27	.972	11	
10	1.00	26545	1.23	25769	1.25	16.4	27	.971	13	
11	1.00	39883	1.00	38761	1.02	16.4	27	.972	14	
12	1.00	33284	1.10	32354	1.11	16.3	27	.972	15	
13	1.00	38947	1.01	37927	1.03	16.6	27	.974	17	
14	1.00	35602	1.06	34699	1.07	16.3	27	.975	18	
15	1.00	17644	1.51	17141	1.53	16.3	27	.971	19	
16	1.00	20121	1.41	19548	1.43	16.3	27	.972	20	
17	1.00	21141	1.38	20512	1.40	16.4	27	.970	22	
18	1.00	3135	3.57	2983	3.66	15.8	27	.952	23	
19	1.00	3549	3.36	3377	3.44	15.9	27	.952	24	
20	1.00	2101	4.36	1982	4.49	15.5	27	.943	25	

100 µl
 10 µl and
 2 µl
 ed on 1

HUE S/F media for 6/23 column Hys Seph + C-18
 1st fraction test 500 µl (out of 1.5 ml fraction) Tube
 dried down resuspended 10 µl - 1 sample
 resuspended 500 µl Tube in 20 µl - 10 µl, 5 µl, 2 µl

HUE S/F media for 6/23 column Hys Seph + C-18
 2nd fraction 500 µl - resuspended in 10 µl - 1 sample

Refer to worksheet of 6/23 for plots of above samples

Mits Asay

1 AFL10 6/24 500ul	2 AFL10 6/24 500ul 20-12	3 AFL10 6/24 500ul 20-6	4 AFL10 6/24 20-9	5 AFL10 6/21 1-fraction 2, 3, 4 12-10ul	6 AFL10 6/21 1-fraction 2, 3, 4 12-2ul
7 HUV S/F HFC dialyzed 500ul	8 HUV S/F HFC dialyzed 500ul	9 HUV S/F C-18 6/23 500ul 12-fac	10 HUV S/F C-18 6/23 500ul 20-12ul	11 HUV S/F C-18 6/23 500ul 20-6ul	12 HUV S/F C-18 6/23 500ul 20-2ul
13 HUV S/F C-18 250ul	14 Cuv 10-4	15 Cuv 5 ng 5 ng	16 Cuv 5 ng	17 Cuv 2 ng	18 Blank
19 Blank	20 Blank	21	22	23	24

G/F on at 5:00 p

Integrase Assay NRK cells (1.4^{plates} / 1.1.1.)

HUVE S/F HAc 10 2nd fraction 6/24 500 μ l Tube
resuspended in 10 μ l 5M HCl - this used as 1 sample
to 500 μ l Tube resuspended in 20 μ l - used as
12 μ l, 5 μ l, 2 μ l, ~~2 μ l~~ samples

HUVE S/F HAc 10 from '1' fraction 2, 3, 4 - 100 μ l
from each fraction combined and dry-saved.
resuspended this in 12 μ l - used as 12 μ l and
2 μ l samples.

HUVE S/F HAc dialyzed 500 μ l dialyzed
media resuspended in 20 μ l 5M HCl - used as 1
sample.

HUVE S/F media from 6/23 column Hep Syn + C-18
1st fraction post 500 μ l (out of 1.5 ml fraction) Tube
dried down resuspended 10 μ l - 1 sample
resuspended 500 μ l Tube in 20 μ l - 12 μ l, 5 μ l, 2 μ l


HUVE S/F media from 6/23 column Hep Syn + C-18
2nd fraction 500 μ l - resuspended in 10 μ l - 1 sample

Refs to western of - for blots of above
samples

KS- DB60 R 32

clones 1, 2, 3, 4



Eco RI R_x cut of
insert



1	1	2	2	new	3	3	4	4
cut	cut	cut	cut	lig	insert	cut	insert	cut
3rd	5th	3rd	5th		3rd	5th	3rd	5th

I picked this clone to grow up + subclone into M13

Southern
Probed w/
λgt11 DB60
fragment



gel

Nitrocellulose
Blot


DB60 R32 - clone selected from Hure & DNH

λ gt11 library of screening w/ DB60
EcoRI fragment picked in 1st screen
w/ 2 P126F antibody

insert cut from λ gt11, ~~for~~ subcloned
into KS+ bluescript phagemid

grown up in 50 ml of NM522 cells,
plasmid prep of 50 ml, ran 10 μ g
on low melt gel + cut out EcoRI
fragment. Ran over Clutip column
(very hard to push low salt wash through,
too much agarose, need smaller slice)

Fragment run on gel w/ markers



Fragment size is